

## Recommended devices for PERIpath Mobile application usage

### iOS\*

#### + PHONES

- Apple iPhone 13
- Apple iPhone SE (3rd gen)

#### + TABLETS

- Apple iPad 2021 (9th gen)
- Apple iPad Pro 11

### Android\*

#### + PHONES

- Samsung Galaxy S21
- Samsung Galaxy S22

#### + TABLETS

- Samsung Galaxy Tab S7

\*Latest operating system version required.

## Minimum specifications for PERIpath Mobile compatible devices

### iOS

|                            |                                |
|----------------------------|--------------------------------|
| System:                    | Latest version                 |
| Device capabilities:       | must be equipped with a camera |
| RAM:                       | 2 GB                           |
| Minimum screen resolution: | 1334 by 750 px                 |

For satisfactory performance, use devices with the A11 SoC or better.

### Android

|                            |   |
|----------------------------|---|
| System:                    | Android 10 or higher (Android 9 is supported but not recommended) |
| Device capabilities:       | must be equipped with a camera                                    |
| RAM:                       | 4 GB  |
| Minimum screen resolution: | 1334 by 750 px  |

For satisfactory performance use at least 2 years old mid-range or 3 years old flagship Android devices.