

Recommended Devices for Using the PERIpath Mobile Application

iOS*

+PHONES

- Apple iPhone 13
- Apple iPhone 15
- Apple iPhone SE (3rd gen)

+TABLETS

- Apple iPad (10th gen)
- Apple iPad Pro 11

Android*

+PHONES

- Samsung Galaxy S21
- Samsung Galaxy S22
- Samsung Galaxy S23

+TABLETS

- Samsung Galaxy Tab S7
- Samsung Galaxy Tab A8
- Samsung Galaxy Tab S8

*Latest operating system version required.

Minimum Specifications for PERIpath Mobile Compatible Devices

iOS

- +System: Latest version of the operating system
- +Device Capabilities: Must be equipped with a camera
- +RAM: Minimum of 3 GB
- +Minimum Screen Resolution: 1334 x 750 pixels

For satisfactory performance, it is recommended to use devices equipped with the A13 SoC or higher.

Android

- +System: Android Version: 12 or higher (Note: The app can be also launched on Android versions 9 to 11, but these versions are no longer supported)
- +Device Capabilities: Must be equipped with a camera
- +RAM: Minimum of 4 GB
- +Minimum Screen Resolution: 1334 x 750 pixels

For satisfactory performance, it is advisable to use Android devices that are at least 2 years old in the mid-range category or 3 years old in the flagship category.